



B3-FO-01-012	Business Development – Engineering	 Technology
	Formulary	
Version: 02	Instruction for Use	

# Implant Tool




English

## Instructions for Use




HL Technology S.A., rue Jardinière 153  
 CH 2300 La Chaux-de-Fonds  
 Tel: +41 (0)32 925 90 50  
[www.hl-technology.ch](http://www.hl-technology.ch)

<b>B3-FO-01-012</b>	<b>Business Development – Engineering</b>	
	<b>Formulary</b>	
<b>Version: 02</b>	<b>Instruction for Use</b>	

## Table of Contents

<b>1</b>	<b>System Description.....</b>	<b>3</b>
1.1	Intended Use .....	3
1.2	Intended user.....	3
1.3	Indications for use and patient population .....	3
1.4	Contra-indications .....	3
1.5	Warnings/ Precautions .....	3
<b>2</b>	<b>Components:.....</b>	<b>4</b>
2.1	Tools.....	4
2.2	Procedure for use .....	5
2.2.1	Screwdriver Adapter Insertion For Manual Use .....	5
2.2.2	Screwdriver Insertion into the Torque Limiter .....	5
<b>3</b>	<b>Packaging / Storage .....</b>	<b>6</b>
<b>4</b>	<b>Treatment:.....</b>	<b>7</b>
4.1	Disinfection.....	7
4.2	Sterilization.....	8
<b>5</b>	<b>Recommendations for use .....</b>	<b>9</b>
<b>6</b>	<b>Disposal .....</b>	<b>9</b>

<b>B3-FO-01-012</b>	<b>Business Development – Engineering</b>	 Technology
	<b>Formulary</b>	
<b>Version: 02</b>	<b>Instruction for Use</b>	

## 1 System Description

Implant tools are intended to tighten and loosen prosthetic screws in laboratories or Dental offices. These screws are used to maintain the abutments in the implants. Temporary or definitive prostheses are fixed on the abutments.

Implant tools or drivers have an ISO or square shanks to be used with an adapted torque wrench to tighten the screws with good torque according to the manufacturer's recommendation.

Implant tools are delivered in kit or individually non-sterile.

### 1.1 Intended Use

Implant Tool is intended to be used with screws and torque wrench to fix and remove dental prosthesis in the mouth.

### 1.2 Intended user

The assembly of the components with the prosthesis must be performed by a properly trained professional.

Cleaning and maintenance of the prosthesis must be advised by the dental practitioner.

### 1.3 Indications for use and patient population

- Implant tools are intended for dental professionals.
- Implant tools are used with screws and torque wrenches.
- Implant tools are used for the replacement or modification of a dental prosthesis or a set of prosthetic elements.
- For fully or partially edentulous jaws.

### 1.4 Contra-indications

Do not use the devices if the drivers or the screws are deteriorated.

### 1.5 Warnings/ Precautions








It is important to check the following points before use.

- Make sure the retreatment (disinfection and sterilization) has been carried out before each use.
- To have chosen the good driver according to the screws.
- To respect the recommendations of the Implant manufacturer concerning the tightening torque.

## 2 Components:

All the components described below can be found as a kit (Implant Buddy) or separately.

### 2.1 Tools

1	Body	
2	Manual Adapter	
3	Hex Screwdriver 0.9 mm	
4	Hex Screwdriver 1.2 mm	
5	Hex Screwdriver 1.25 mm	
6	Hex Screwdriver 1.3 mm	
7	Angled Screwdriver	
8	SCS Screwdriver	
9	Unigrip Screwdriver	











FIG.1 : Implant Buddy Dental

1	Body	
2	Manual Adapter	
3	Multi-Unit Driver	
4	Hex Screwdriver 0.9 mm	
5	Hex Screwdriver 1.2 mm	
6	Hex Screwdriver 1.3 mm	
7	Angled Screwdriver	
8	Straight Driver	
9	SCS Screwdriver	
10	Unigrip Screwdriver	





FIG.2 : Implant Buddy Latch

B3-FO-01-012	Business Development – Engineering	 Technology
	Formulary	
Version: 02	Instruction for Use	

## 2.2 Procedure for use

During tightening operations, the screws must be tightened with a torque limiter that is compatible with the screwdriver.

For all use of the device, its components must be used according to the method described below:

### 2.2.1 Screwdriver Adapter Insertion For Manual Use




FIG.3 : Screwdriver Adapter Mounting



FIG.4 : Screwdriver Adapter Mounted




FIG.5 : Screwdriver Adapter Mounting




FIG.6 : Screwdriver Adapter Mounted

### 2.2.2 Screwdriver Insertion into the Torque Limiter





FIG.7: Mounting the screwdriver into the torque limiter




FIG.8: The screwdriver in the mounted position

B3-FO-01-012	Business Development – Engineering	 Technology
	Formulary	
Version: 02	Instruction for Use	

### 3 Packaging / Storage

The parts are delivered in a vacuum sealed plastic bag. If the pouch is no longer sealed or under vacuum upon receipt of the parts, they must be returned to the distributor. The kits Implant Buddy are delivered in tube / rigid plastic packaging.

The parts should be kept in a clean, dry place and protected from direct sunlight. The temperature of the storage conditions must remain at room temperature.





Left picture: plastic bag for drivers. Right picture: tube for supports



## 4 Treatment:

### 4.1 Disinfection




<b>B3-FO-01-012</b>	<b>Business Development – Engineering</b>	 Technology
	<b>Formulary</b>	
<b>Version: 02</b>	<b>Instruction for Use</b>	

## 4.2 Sterilization

<b>US Market</b>	
Method	Steam sterilization according to the EN ISO 17665 norm and ANSI/AAMI ST79
Cycle	
1) Pre-vacuum	
2) Sterilisation	4 minutes à 132°C / 270°F
3) Drying time	60 minutes
4) Cooling time	60 minutes at ambient temperature

<b>EU Market</b>	
Method	Steam sterilization according to the EN ISO 17665 norm
Cycle	
1) Pre-vacuum	
2) Sterilisation	18 minutes at 134°C / 273°F at 2 bars
3) Drying time	20 minutes



<b>B3-FO-01-012</b>	<b>Business Development – Engineering</b>	 Technology
	<b>Formulary</b>	
<b>Version: 02</b>	<b>Instruction for Use</b>	

## 5 Recommendations for use

- It is recommended to use the screwdrivers always with a torque wrench to avoid over or under tightening the abutments screws.
- Always following the tightening recommendations of the implant manufacturers.
- Never exceed the tightening values in the table below.

<b>Driver</b>	<b>Torque max</b>
Hex Screwdriver 0.9 mm	<b>15N/cm</b>
Hex Screwdriver 1.2 mm	<b>32N/cm</b>
Hex Screwdriver 1.25 mm	<b>30N/cm</b>
Hex Screwdriver 1.3 mm	<b>30N/cm</b>
SCS Screwdriver	<b>35N/cm</b>
Unigrip Screwdriver	<b>35N/cm</b>
Straight Driver	<b>15N/cm</b>
Angled Driver	<b>30N/cm</b>

## 6 Disposal

Disposal must be done in accordance with the regulations applied in the country of use.